

## Product Overview

- **Articulating Block**  
Allows the cage to safely and accurately pivot to final position.


- **Pyramidal teeth**  
Enhances implant stability.

- **Large graft cavity**  
Accommodates bone graft packing to aid in the fusion process.


- **Self-distracting nose**  
Allows for easier of insertion.

- **Radiographic markers**  
Produces clear radiographic visualization during implant insertion and in final position.

- **Anterior and posterior side holes**  
Contributes to the vascularization of the graft.




- **Lordotic angle**  
Provides extensive surface contact and restores the spine's natural lordotic curve.




The Velofix™ TLIF Cage consists of implants available in various heights, widths, lengths and lordotic angles with an open architecture to accept packing of bone graft material. The device consists of the cage body, the articulating block, and radiographic markers. The articulating block allows the cage to pivot to the final positioning.

The Velofix™ TLIF Cage may be implanted as a single device via a transforaminal approach, and can be adjusted to align with the anterior anatomy of the vertebral endplates.



### PATIENT POSITIONING AND SURGICAL EXPOSURE

Identify the affected disc with fluoroscopy. Expose the intervertebral space using the surgeon's customary transforaminal approach.




**Fig 1**

### FACET RESECTION

In order to gain access to the disc, resect the facet joints on the desired side for implant insertion(**Fig 2**).

*Note: The facet joint can be preserved if sufficient access exists for discectomy and implant insertion or if the approach is extraforaminal.*



**Fig 2**

### Pedicle screw installation

Install a pedicle screw systems (Perfix™ / Perfix™ MIS, ANAX™ 5.5 / ANAX™ 5.5 MIS) in the levels above and below the affected disc. Screws can be used for distraction, which may be necessary for the discectomy, as well as the trial and implantation steps.




Fig 3

### DISTRACTION OF THE DISC SPACE

Attach the T-Handle to the Distractor.

Insert Distractor, then rotate to restore the natural disc height(Fig 4). The distractors are available in heights from 7-16mm in 1mm increments.




Fig 4

|               | Instrument           |
|---------------|----------------------|
| TL0380        | Distractor 7mm       |
| CB0108-CB0116 | S. Distractor 8-16mm |
| SP0002        | T-Handle             |

*Note: In the case of an extremely collapsed disc, a rod can be installed on the side opposite implant insertion. Once preliminarily locked, the rod will help maintain distraction.*

## DISC SPACE PREPARATION

Connect the T-Handle with the appropriate Bone Shaver. Insert into the vertebral body and rotate to remove residual intradiscal material and to create a channel in the dorsal most endplate(**Fig 5**).

### DISC SHAVERS


|               | Instrument         |
|---------------|--------------------|
| PL0310-PL0390 | Bone Shaver 7-16mm |
| SP0002        | T-Handle           |

### CURETTES and BONE RASP


|        | Instrument                  |
|--------|-----------------------------|
| TL0400 | Ring Curette – Straight     |
| TL0410 | Ring Curette – Right Curved |
| TL0420 | Ring Curette – Left Curved  |
| TL0440 | Bone Rasp                   |

Remove the remaining soft tissue and cartilaginous endplate coverings using the Curettes and Bone Rasp.

**(Fig 6)**



**Fig 5**



**Fig 6**

## ASSEMBLY THE TRIAL

### 1. TRIAL ASSEMBLY

Attach the Rotating Holder to the Trial Inserter with the correct orientation(**Fig 7**).

Select the appropriate trial head for the disc size and attach to the Trial Inserter(**Fig 8**).

To lock, press the button on the Rotating Holder and push so the Trial Head makes contact with the Rotating Holder.

*Note: The Trial Heads have the same heights as the implants and make it possible to determine the appropriate height.*




Fig 7




Fig 8

|               | Instrument        |
|---------------|-------------------|
| TL0060        | Trial Inserter    |
| TL0070-TL0160 | Trial Head 7-16mm |
| TL0230        | Rotating Holder   |

### 2. SIZING THE DISC SPACE

Lock the Trial Head on the Trial Inserter and insert the Trial Head into the disc space until it contacts the fibrous annulus. (**Fig 9**)





Fig 9

|               | Instrument        |
|---------------|-------------------|
| TL0060        | Trial Inserter    |
| TL0070-TL0160 | Trial Head 7-16mm |
| TL0230        | Rotating Holder   |

**3. FINAL POSITIONING**


Unlock the Rotating Holder(Fig 11) and advance the trial to final position(Fig 10). Confirm the position by fluoroscopy.

Remove the Trial assembly using the Sliding Hammer. (Fig 12)




**Fig 10**

|        | Instrument     |
|--------|----------------|
| PL0430 | Sliding Hammer |



**Fig 11 Locking and Unlocking the Rotating Holder**



**Fig 12 Assembling the Sliding Hammer**

### ATTACH CAGE TO INSERTER

Attach the Rotating Holder to the inserter. Set it to the LOCK position by pressing the Rotating Holder button.

Attach the appropriate Velofix™ TLIF Cage to the inserter(**Fig 13**). Insert the Fixation Bolt into the Inserter and rotate to tighten(**Fig 14**).

Place the Rotating Holder in the locked position by pressing the button on the Rotating Holder.




Fig 13

|        | Instrument             |
|--------|------------------------|
| TL0210 | Fixation Bolt-Straight |
| TL0220 | Inserter-Straight      |
| TL0230 | <u>Rotating Holder</u> |

\* Note : Before assembly the Inserter with Cage, Check the direction of the Articulating Block.(**Fig 13**)




Fig 14

### BONE PACKING

Place the Velofix™ TLIF Cage on the Packing Platform and carefully compact graft using a Bone Packing Bar (**Fig 15**). Fill until graft is level with the top of the cage to ensure optimal contact with the vertebral endplates.

*Note: Autologous bone is options for filling the Velofix™ TLIF Cage.*




Fig 15

|        | Instrument       |
|--------|------------------|
| TL0520 | Packing Platform |
| PL0450 | Bone Packing Bar |

### IMPLANT INSERTION

Insert the Velofix™ TLIF Cage into the disc space until it contacts the fibrous annulus(**Fig 16**).

Press the button on the Rotating Holder and gently pull upwards to unlock(**Fig 17**), and advance the implant to final position.

Confirm the position by fluoroscopy.

|        | Instrument             |
|--------|------------------------|
| TL0210 | Fixation Bolt-Straight |
| TL0220 | Inserter-Straight      |
| TL0230 | <u>Rotating Holder</u> |




Fig 16




Fig 17

### INSTALL POSTERIOR GRAFT

After Velofix™ TLIF Cage is implanted, fill the posterior disc space and the lateral disc space with bone graft for fusion using the Impactor(**Fig 18**).

|        | Instrument |
|--------|------------|
| TL0510 | Impactor   |




Fig 18



**STABILIZATION AND COMPRESSION**

Once the Velofix™ TLIF Cage is in place, compress the segment with posterior screws and rod to stabilize the Velofix™ TLIF Cage(Fig 19).




**Fig 19**

**IMPLANT REMOVAL**

Should removal or revision be deemed necessary by the surgeon, the implant can be removed by using the Inserter and the Fixation Bolt. Attach the Inserter and Fixation Bolt to the Velofix™ TLIF Cage and remove the cage frame from the disc space. If greater force is needed, attach the Sliding Hammer to the Fixation Bolt. Gently impact the Sliding Hammer to remove the implant(Fig 20).

|        | Instrument             |
|--------|------------------------|
| TL0210 | Fixation Bolt-Straight |
| TL0220 | Inserter-Straight      |
| PL0430 | Sliding Hammer         |



**Fig 20**

## Order information

### Implants (Single Use Only)



| No.                    | Part No. | Part Description              |
|------------------------|----------|-------------------------------|
| <b>TLIF Cage L28mm</b> |          |                               |
| 1                      | TLA1007  | TLIF CAGE L28mm X A0° X H7mm  |
| 2                      | TLA1008  | TLIF CAGE L28mm X A0° X H8mm  |
| 3                      | TLA1009  | TLIF CAGE L28mm X A0° X H9mm  |
| 4                      | TLA1010  | TLIF CAGE L28mm X A0° X H10mm |
| 5                      | TLA1011  | TLIF CAGE L28mm X A0° X H11mm |
| 6                      | TLA1012  | TLIF CAGE L28mm X A0° X H12mm |
| 7                      | TLA1013  | TLIF CAGE L28mm X A0° X H13mm |
| 8                      | TLA1014  | TLIF CAGE L28mm X A0° X H14mm |
| 9                      | TLA1015  | TLIF CAGE L28mm X A0° X H15mm |
| 10                     | TLA1016  | TLIF CAGE L28mm X A0° X H16mm |
| 11                     | TLA1507  | TLIF CAGE L28mm X A5° X H7mm  |
| 12                     | TLA1508  | TLIF CAGE L28mm X A5° X H8mm  |
| 13                     | TLA1509  | TLIF CAGE L28mm X A5° X H9mm  |
| 14                     | TLA1510  | TLIF CAGE L28mm X A5° X H10mm |
| 15                     | TLA1511  | TLIF CAGE L28mm X A5° X H11mm |
| 16                     | TLA1512  | TLIF CAGE L28mm X A5° X H12mm |
| 17                     | TLA1513  | TLIF CAGE L28mm X A5° X H13mm |
| 18                     | TLA1514  | TLIF CAGE L28mm X A5° X H14mm |
| 19                     | TLA1515  | TLIF CAGE L28mm X A5° X H15mm |
| 20                     | TLA1516  | TLIF CAGE L28mm X A5° X H16mm |
| 21                     | TLA1807  | TLIF CAGE L28mm X A8° X H7mm  |
| 22                     | TLA1808  | TLIF CAGE L28mm X A8° X H8mm  |
| 23                     | TLA1809  | TLIF CAGE L28mm X A8° X H9mm  |
| 24                     | TLA1810  | TLIF CAGE L28mm X A8° X H10mm |
| 25                     | TLA1811  | TLIF CAGE L28mm X A8° X H11mm |
| 26                     | TLA1812  | TLIF CAGE L28mm X A8° X H12mm |
| 27                     | TLA1813  | TLIF CAGE L28mm X A8° X H13mm |
| 28                     | TLA1814  | TLIF CAGE L28mm X A8° X H14mm |

|                        |         |                               |
|------------------------|---------|-------------------------------|
| 29                     | TLA1815 | TLIF CAGE L28mm X A8° X H15mm |
| 30                     | TLA1816 | TLIF CAGE L28mm X A8° X H16mm |
| <b>TLIF Cage L32mm</b> |         |                               |
| 31                     | TLA2007 | TLIF CAGE L32mm X A0° X H7mm  |
| 32                     | TLA2008 | TLIF CAGE L32mm X A0° X H8mm  |
| 33                     | TLA2009 | TLIF CAGE L32mm X A0° X H9mm  |
| 34                     | TLA2010 | TLIF CAGE L32mm X A0° X H10mm |
| 35                     | TLA2011 | TLIF CAGE L32mm X A0° X H11mm |
| 36                     | TLA2012 | TLIF CAGE L32mm X A0° X H12mm |
| 37                     | TLA2013 | TLIF CAGE L32mm X A0° X H13mm |
| 38                     | TLA2014 | TLIF CAGE L32mm X A0° X H14mm |
| 39                     | TLA2015 | TLIF CAGE L32mm X A0° X H15mm |
| 40                     | TLA2016 | TLIF CAGE L32mm X A0° X H16mm |
| 41                     | TLA2507 | TLIF CAGE L32mm X A5° X H7mm  |
| 42                     | TLA2508 | TLIF CAGE L32mm X A5° X H8mm  |
| 43                     | TLA2509 | TLIF CAGE L32mm X A5° X H9mm  |
| 44                     | TLA2510 | TLIF CAGE L32mm X A5° X H10mm |
| 45                     | TLA2511 | TLIF CAGE L32mm X A5° X H11mm |
| 46                     | TLA2512 | TLIF CAGE L32mm X A5° X H12mm |
| 47                     | TLA2513 | TLIF CAGE L32mm X A5° X H13mm |
| 48                     | TLA2514 | TLIF CAGE L32mm X A5° X H14mm |
| 49                     | TLA2515 | TLIF CAGE L32mm X A5° X H15mm |
| 50                     | TLA2516 | TLIF CAGE L32mm X A5° X H16mm |
| 51                     | TLA2807 | TLIF CAGE L32mm X A8° X H7mm  |
| 52                     | TLA2808 | TLIF CAGE L32mm X A8° X H8mm  |
| 53                     | TLA2809 | TLIF CAGE L32mm X A8° X H9mm  |
| 56                     | TLA2810 | TLIF CAGE L32mm X A8° X H10mm |
| 55                     | TLA2811 | TLIF CAGE L32mm X A8° X H11mm |
| 56                     | TLA2812 | TLIF CAGE L32mm X A8° X H12mm |
| 57                     | TLA2813 | TLIF CAGE L32mm X A8° X H13mm |
| 58                     | TLA2814 | TLIF CAGE L32mm X A8° X H14mm |
| 59                     | TLA2815 | TLIF CAGE L32mm X A8° X H15mm |
| 60                     | TLA2816 | TLIF CAGE L32mm X A8° X H16mm |

### Instruments

TL0210 Fixation Bolt-Straight



TL0220 Inserter-Straight



TL0230 Rotating Holder



TL0060 Trial Inserter



TL0070 Trial Head 7mm



TL0080 Trial Head 8mm



TL0090 Trial Head 9mm



TL0100 Trial Head 10mm



TL0110 Trial Head 11mm



TL0120 Trial Head 12mm




TL0130 Trial Head 13mm



TL0140 Trial Head 14mm



TL0150 Trial Head 15mm



TL0160 Trial Head 16mm



TL0400 Ring Curette – Straight



TL0410 Ring Curette – Left Curved



TL0420 Ring Curette – Right Curved




TL0440 Bone Rasp



PL0450 Bone Packing Bar



TL0520 Packing Platform



TL0510 Impactor



PL0310 Bone Shaver - 7mm



PL0320 Bone Shaver - 8mm



PL0330 Bone Shaver - 9mm



PL0340 Bone Shaver - 10mm



PL0360 Bone Shaver - 12mm



PL0380 Bone Shaver - 14mm



PL0390 Bone Shaver - 16mm



TL0380 Distractor - 7mm



CB0108 S. Distractor - 8mm



CB0109 S. Distractor - 9mm



CB0110 S. Distractor - 10mm



CB0112 S. Distractor - 12mm



CB0114 S. Distractor - 14mm




CB0116 S. Distractor - 16mm



PL0430 Sliding Hammer



SP0002 T-handle



**◆ Intended Use / Indications For Use**

The Velofix™ TLIF Cage is indicated for intervertebral body fusion at one level or two contiguous levels in the lumbar spine from L2 to S1 in patients with degenerative disc disease (DDD) with up to Grade I spondylolisthesis at the involved level(s). DDD is defined as back pain of discogenic origin with degeneration of the disc confirmed by patient history and radiographic studies. These patients should be skeletally mature and have had six months of non-operative treatment.

The Velofix™ TLIF Cage is to be used with supplemental internal spinal fixation. Additionally, the Velofix™ TLIF Cage is to be used with autogenous bone graft material.

**◆ Contraindications**

Contraindications include, but are not limited to:

- The Velofix™ TLIF Cage should not be implanted in patients with an active infection at the operative site.
- Any active or suspected latent infection in or about the spine.
- Any mental or neuromuscular disorder which would create an unacceptable risk of fixation failure or complications in postoperative care.
- Bone stock compromised by disease, infection or prior implantation which cannot provide adequate support and/or fusion to the devices.
- Obesity. An overweight or obese patient can produce loads on the spinal system which can lead to failure of the fusion of the device or to failure of the device itself. Obesity is defined according to the W.H.O. standards.
- Recent infection, fever, or leukocytosis.
- Open wounds.
- Metal and / or foreign body sensitivity, documented or suspected.
- Patient having inadequate tissue coverage over the operative site.
- Pregnancy
- Excessive local inflammation.
- Prior fusion at the level(s) to be treated.
- Other medical or surgical conditions which would preclude the potential benefit of spinal implant surgery, such as the presence of tumors, congenital abnormalities, elevation of sedimentation rate unexplained by other diseases, elevation of white blood cell count(WBC), or marked left shift in the WBC differential count.
- A condition of senility, mental illness, or substance abuse. These conditions, among others, may cause the patient to ignore certain necessary limitations and precautions in the use of the implant, leading to failure or other complications.
- Any neuromuscular deficit which places an unsafe load level on the device during the healing period.
- Any abnormality present which affects the normal process of bone remodeling including, but not limited to,

severe osteoporosis involving the spine, bone absorption, osteopenia, primary or metastatic tumors involving the spine, active infection at the site or certain metabolic disorders affecting osteogenesis.

#### ◆ ⚠ Warnings

- Benefit of spinal fusions utilizing any interbody fusion system has not been adequately established in patients with stable spines.
- Potential risks associated with the use of this system, which may require additional surgery, include: device component fracture, loss of fusion, non-union, fracture of the vertebra, neurological injury, and vascular or visceral injury.
- Discard all damaged or mishandled implants.
- Never reuse an implant, even though it may appear undamaged. Reusing of the implant can cause early fracture, deformation and patient's infection.
- Interbody fusion devices cannot withstand activity and load levels equal to those placed on normal healthy bone. Until maturation of the fusion mass is confirmed, do not subject this device to the stress of full weight bearing, or implant failure may result.
- Mixing Metal: Some degree of corrosion occurs on all implanted metal and alloys. Contact of dissimilar metals, however, may accelerate this corrosion process. The presence of corrosion may accelerate fatigue fracture of implants, and the amount of metal compounds released into the body system will also increase.
- Because different manufacturers employ different materials, varying tolerances and manufacturing specifications, and differing design parameters, the Velofix™ TLIF Cage should not be used in conjunction with components from any other manufacturer's spinal system and instruments. Any such use will negate the responsibility of U&I Corporation for the performance of the resulting mixed component implant.
- Any decision by a physician to remove the interbody fusion device should take into consideration such factors as the risk to the patient of the additional surgical procedure as well as the difficulty of removal.
- Implant removal should be followed by adequate postoperative management to avoid fracture.

#### ◆ Precautions

**PHYSICIAN NOTE:** Although the physician is the learned intermediary between the company and the patient, the important medical information given in this document should be conveyed to the patient.

**CAUTION: FEDERAL LAW (USA) RESTRICTS THESE DEVICES TO SALE BY OR ON THE ORDER OF A PHYSICIAN**

- Based on the fatigue testing results, the physician/surgeon should consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc. which may have impact on the performance of the system.

- Patients who smoke have been shown to have an increased incidence of non-unions. Such patients should be advised of this fact and warned of the potential consequences.
- If the patient is involved in an occupation or activity which applies inordinate stress upon the implant (e.g., substantial walking, running, lifting, or muscle strain) resultant forces can cause failure of the device.
- Safety and effectiveness have not been established for patients with the following conditions: previous fusion attempt at the involved level(s), spondylolisthesis greater than Grade I, three or more levels to be fused, concomitant conditions requiring steroids, systemic or terminal illness, active drug abuse, gross obesity and severe osteoporotic conditions.
- In some cases, the progression of degenerative disease may be so advanced at the time of implantation that it may substantially decrease the expected useful life of the appliance. For such cases, orthopedic devices can only be considered a delaying technique or temporary remedy.
- Before clinical use, the surgeon should thoroughly understand all aspects of the surgical procedure and limitations of the interbody fusion. Requiring detailed knowledge of spinal surgery. This device is recommended for use only by surgeons familiar with preoperative and surgical techniques, cautions, and potential risks associated with such spinal surgery. Knowledge of surgical techniques, proper reduction, selection and placement of implants, and pre-/post-operative patient management are considerations essential to a successful surgical outcome.
- Patients should be instructed in detail about the limitations of the implants, including, but not limited to, the impact of excessive loading through patient weight or activity, and be taught to govern their activities accordingly. The patient should understand that a metallic or polymer implant is not as strong as normal, healthy bone and will bend, loosen or fracture if excessive demands are placed on it. An active, debilitated, or demented patient who cannot properly use weight supporting devices may be particularly at risk during postoperative rehabilitation.
- Appropriate selection, placement and fusion of the spinal system components are critical factors which affect implant service life. Also, as in the case of all prosthetic implants, the durability of these components is affected by numerous biological, biomechanical and other extrinsic factors, which limit their service life. Accordingly, strict adherence to the indications, contraindications, precautions, and warnings for this product is essential to potentially maximize service life.

***Note: While proper implant selection can minimize risks, the size and shape of human bones present limitations on the size, shape, and strength of the implants.***

- Care must be taken to protect the components from being marred, nicked or notched as a result of contact with metal or abrasive objects. Alterations will produce defects in surface finish and internal stresses which may become the focal point for eventual breakage of the implant.
- Check if implant type is appropriate for use in intended region of spine.
- Patients with previous spine surgery at the level(s) to be treated may have different clinical outcomes compared to those without a previous surgery.

△ **CAUTION: For use on or by the order of a physician only.**

△ **Federal law(USA) restricts these devices to sale by or on the order of a physician.**